

2017-F3-11 Equipment set for thermometer calibration

Contact person: NAKATANI Masao

Information of Equipment:

Fluke DAQ (1586A) with Scanner, Dry Well (9142), Micro Bath (7103) etc.
<http://us.flukecal.com/products/temperature-calibration/digital-thermometer-readouts>

Temperature range is -30 degC to 150 degC.

You can calibrate some thermometers at the same time.

Condition of Use and Remarks

Must contact with the responsible person prior to the application.


Operate the equipment by yourself in principle.

Application periods:

Any time, as needed

Number of retained
equipment:

1 set



[Contact Information]

1-1-1 Yayoi, Bunkyo-ku, Tokyo 113-0032

Earthquake Research Institute, the University of Tokyo

Research Support Team (Joint Usage Section)

Phone: 03-5841-5710, 1769 FAX: 03-5689-4467

Email: k-kyodoriyo@eri.u-tokyo.ac.jp